

Zak Rabbit Pab - Zak

Item Number: CM0016193

Introduction

Rabbit anti-human ZAK polyclonal antibody validated for Western blot (WB), immunofluorescence (IF/ICC), and ELISA. Cross-reacts with mouse and rat ZAK. Useful for investigating ZAK kinase function in disease.

[Learn More](#)

Parameter	Specification
Product Name	ZAK Rabbit pAb
Remarks/Alias	pk, AZK, MLT, MRK, ZAK, CNM6, MLK7, MLTK, SFMMP, mlklak, MLTKbeta, MLTKalpha
Species	Human
Gene ID (Human)	51776
Gene ID	51776
Immunogen	Recombinant protein; recombinant fusion protein containing a sequence corresponding to amino acids 100-312 of human ZAK (NP_598407.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
Protein Weight	35kDa/51kDa/91kDa
Shipping	Ice bag