


Zap70 Rabbit Monoclonal Antibody For Wb, If-P, Elisa - P43403

Item Number: CM0005184



Introduction

Recombinant rabbit monoclonal antibody against human ZAP70, validated for Western blot (WB), immunofluorescence (IF-P), and ELISA. Target tyrosine kinase essential for T-cell receptor signaling and adaptive immunity. 70 kDa.

[Learn More](#)

Parameter	Value
Product Name	ZAP70 Rabbit mAb
Remarks/Alias	SRK; STD; TZK; STCD; IMD48; ADMIO2; ZAP-70; ZAP70
Species	Human
GeneID (Human)	7535
GeneID	7535
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 550-619 of human ZAP70 (NP_001070.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human
SWISS	P43403
Protein Weight	70kDa
Shipping	Ice bag