


Zeb1 Rabbit Monoclonal Antibody For Wb, Elisa - P37275

Item Number: CM0007320



Introduction

ZEB1 Rabbit Monoclonal Antibody (P37275) is a high-quality tool for WB and ELISA, with reactivity in human and rat samples. Essential for studying EMT, transcriptional repression, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	ZEB1 Rabbit mAb
Remarks/Alias	BZP; TCF8; AREB6; FECD6; NIL2A; PPCD3; ZFHEP; ZFHX1A; DELTAEF1; ZEB1
Species	Human
Gene ID (Human)	6935
Gene ID	6935
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 800-900 of human ZEB1 (NP_110378.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P37275
Protein Weight	124kDa
Shipping	Ice bag