

Znf131 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P52739

Item Number: CM0015899

Introduction

Rabbit polyclonal ZNF131 antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, and rat. ZNF131 is a zinc finger transcriptional repressor involved in brain function and development.

[Learn More](#)

Parameter	Value
Product Name	ZNF131 Rabbit pAb
Remarks/Alias	ZBTB35; pHZ-10; ZNF131
Species	Human
GeneID (Human)	7690
GeneID	7690
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-240 of human ZNF131 (NP_003423.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P52739
Protein Weight	71kDa
Shipping	Ice bag