

Znf148 Rabbit Pab - Q9Uqr1

Item Number: CM0026212

Introduction

ZNF148 Rabbit pAb is a polyclonal antibody raised against recombinant ZNF148 protein. It is suitable for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for studying transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	ZNF148 Rabbit pAb
Remarks/Alias	BERF-1; BFCOL1; GDACCF; ZBP-89; ZFP148; pHZ-52; HT-BETA; ZNF148
Species	Human
GeneID (Human)	7707
GeneID	7707
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 690-794 of human ZNF148 (NP_068799.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQR1
Protein Weight	89kDa
Shipping	Ice bag