

Znf24 Rabbit Polymab® Monoclonal Antibody For Wb, If/Icc, If-P, Ihc-P, Elisa - P17028

Item Number: CM0009200



Introduction

Rabbit recombinant monoclonal antibody to human ZNF24 (KOX17, ZNF191, ZSCAN3, RSG-A) for WB, IF/ICC, IF-P, IHC-P, ELISA; cross-reacts with mouse, rat. ~42kDa. Ideal for transcription factor and myelination studies.

[Learn More](#)

Parameter	Value
Product Name	ZNF24 Rabbit PolymAb®
Remarks/Alias	KOX17; RSG-A; ZNF191; ZSCAN3; Zfp191; ZNF24
Species	Human
Gene ID (Human)	7572
Gene ID	7572
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 160-242 of human ZNF24 (NP_110378.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17028
Protein Weight	42kDa
Shipping	Ice bag