

Znf384 Rabbit Pab For Wb, Ihc-P, Elisa - Q8Tf68

Item Number: CM0023879

Introduction

CamelBio offers IVD-grade ZNF384 Rabbit Polyclonal Antibody validated for WB, IHC-P, ELISA in human, mouse, rat. Nucleus-localized transcription factor, 63kDa, UniProt Q8TF68.

[Learn More](#)

Parameter	Value
Product Name	ZNF384 Rabbit pAb
Remarks/Alias	NP; CIZ; NMP4; CAGH1; ERDA2; TNRC1; CAGH1A; ZNF384
Species	Human
GeneID (Human)	171017
GeneID	171017
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8TF68
Protein Weight	63kDa
Shipping	Ice bag