

Znf444 Rabbit Pab - Q8N0Y2

Item Number: CM0023988

Introduction

ZNF444 Rabbit pAb detects ZNF444 (EZF2) across human, mouse, and rat samples via WB, IHC-P, IF/ICC, ELISA. A transcriptional regulator binding the SCARF1 promoter, with an observed molecular weight of 35 kDa.

[Learn More](#)

Parameter	Value
Product Name	ZNF444 Rabbit pAb
Remarks/Alias	EZF2; EZF-2; ZSCAN17; ZNF444
Species	Human
Gene ID	55311
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8N0Y2
Protein Weight	35kDa
Shipping	Ice bag