

Znf502 Rabbit Pab - Q8Tbz5

Item Number: CM0014996

Introduction

ZNF502 rabbit polyclonal antibody for Western blot (WB) and ELISA. ZNF502 is a zinc finger protein potentially involved in transcriptional regulation. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	ZNF502 Rabbit pAb
Remarks/Alias	ZNF502
Species	Human
GeneID (Human)	91392
GeneID	91392
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human ZNF502 (NP_149987.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8TBZ5
Protein Weight	63kDa
Shipping	Ice bag