

Znf544 Rabbit Pab - Q6Nx49

Item Number: CM0029223

Introduction

ZNF544 Rabbit pAb for WB and ELISA, detects human and mouse ZNF544. This 82 kDa polyclonal antibody targets the N-terminal region, suitable for transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	ZNF544 Rabbit pAb
Remarks/Alias	ZNF544
Species	Human
GeneID (Human)	27300
GeneID	27300
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6NX49
Protein Weight	82kDa
Shipping	Ice bag