


Zo-1 Rabbit Pab For Wb, If/Icc, Elisa - Q07157

Item Number: CM0016711



Introduction

Rabbit polyclonal antibody against human ZO-1 (TJP1). Validated for WB, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for tight junction and epithelial barrier research.

[Learn More](#)

Parameter	Value
Product Name	ZO-1 Rabbit pAb
Remarks/Alias	ZO-1
Species	Human
GeneID (Human)	7082
GeneID	7082
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-273 of human TJP1 (NP_003248.3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q07157
Protein Weight	195kDa
Shipping	Ice bag