

Zpr1 Rabbit Polyclonal Antibody For Wb And Elisa - 075312

Item Number: CM0014849

Introduction

ZPR1 rabbit polyclonal antibody validated for WB and ELISA. Detects human ZPR1, a zinc finger protein involved in cell cycle progression and SMN1 localization. Ideal for neuroscience and cell signaling research.

[Learn More](#)

Parameter	Value
Product Name	ZPR1 Rabbit pAb
Remarks/Alias	GKAF; ZNF259; ZPR1
Species	Human
Gene ID	8882
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 110-459 of human ZPR1 (NP_003895.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O75312
Protein Weight	51kDa
Shipping	Ice bag