

Zscan4 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Nam6

Item Number: CM0018391

Introduction

Rabbit polyclonal antibody targeting human ZSCAN4 (ZNF494), generated against recombinant protein (aa 1-227). Suitable for Western blot and ELISA, ideal for stem cell pluripotency and telomere dynamics studies.

[Learn More](#)

Parameter	Value
Product Name	ZSCAN4 Rabbit pAb
Remarks/Alias	ZNF494; ZSCAN4
Species	Human
Gene ID (Human)	201516
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-227 of human ZSCAN4 (NP_689890.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8NAM6
Protein Weight	49kDa
Shipping	Ice bag