

Zscan4C Rabbit Pab - 245109

Item Number: CM0014281

Introduction

ZSCAN4C Rabbit pAb is a polyclonal antibody targeting mouse ZSCAN4, with reactivity in mouse. Ideal for WB and ELISA, detecting a 57kDa band. A key tool for stem cell and developmental biology research, supporting studies on pluripotency and early embryogenesis.

[Learn More](#)

Parameter	Value
Product Name	ZSCAN4C Rabbit pAb
Remarks/Alias	Gm397, Zscan4d, XM_142517, ZSCAN4C
Species	Mouse
Gene ID (Human)	201516
Gene ID (Mouse)	245109
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 400-506 of mouse ZSCAN4 (NP_001013787.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
Protein Weight	57kDa
Shipping	Ice bag